

## System Save

1. You must have all the tapes that you plan to use in this save initialized and ready to use prior to starting the save procedure.
  - a. To determine how many tapes you will need you should do the following.
    - i. Type WRKSYSSTS<Enter>
    - ii. Look at “System ASP” . . . . . : 12.58 G
    - iii. Look at “% system ASP used” . . : 64.6923
    - iv. Use rounding, take 13 G – 64% =
2. Must have a dedicated system. Everyone off except you.
3. You must do the save from the console terminal.
4. After everyone is off type in OEKILL<Enter>
5. Look in WRKACTJOB and verify OECNTL is not running under QBATCH.

On a command line type:  
GO SAVE<Enter>

Roll up or page down and select option 21. (System Save)

You will get the screen “Specify Command Defaults.  
Change Prompt for commands to “N”.  
Change Check for active files to “N”.  
Change Message Queue delivery to \*NOTIFY  
Leave the rest the same.

If/when you fill a tape you will receive a message to put in next tape and take the “G” option to go.

Your save may take anywhere from 3 to 8 hours to complete depending on the speed of your processor and the amount of data you are saving.